

Using acceleration sensors to identify rigidity release threshold during Deep Brain Stimulation surgery

Ashesh Shah, Jerome Coste, Jean-Jacques Lemaire, Erik Schkommodau, Raphael Guzman, Ethan Taub, Simone Hemm-Ode

► To cite this version:

Ashesh Shah, Jerome Coste, Jean-Jacques Lemaire, Erik Schkommodau, Raphael Guzman, et al.. Using acceleration sensors to identify rigidity release threshold during Deep Brain Stimulation surgery. 7th international IEEE EMBS Conference on Neural Engineering, Apr 2015, Montpellier, France. poster 520, 2015, IEEE/EMBS Conference on Neural Engineering (NER), 2015. hal-01866559

HAL Id: hal-01866559 https://uca.hal.science/hal-01866559

Submitted on 19 Sep 2018 $\,$

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.





USING ACCELERATION SENSORS TO IDENTIFY RIGIDITY RELEASE THRESHOLD **DURING DEEP BRAIN STIMULATION SURGERY**

A. Shah¹, Member, IEEE, J. Coste², JJ. Lemaire², E. Schkommodau¹, Member, IEEE , R. Guzman³, E. Taub³ and S. Hemm-Ode¹, Member, IEEE ¹Institute for Medical and Analytical Technologies, University of Applied Sciences and Arts Northwestern Switzerland, School of Life Sciences, Muttenz, Switzerland ²CHU de Clermont-Ferrand, EA 7282, IGCNC, Université d'Auvergne, France, CHU de Clermont-Ferrand, France ³Departments of Neurosurgery and Biomedicine, University Hospital Basel, Basel, Switzerland

Background

Deep brain stimulation (DBS) is now a widely accepted surgical treatment for Parkinson's disease (PD). Electrodes are implanted in the patient's brain after intraoperative test stimulation. Changes in parkinsonian rigidity during test stimulation are detected by an evaluator, usually a neurologist, by identifying changes in the resistance of the patient's arm to a passive movement. When a stimulation-induced reduction in rigidity is observed, the stimulation amplitude is noted; this is the clinical rigidity release threshold. The aim of the present study was to test the hypothesis that, at the moment of reduction in rigidity, the speed with which the evaluator moves the patient's arms increases, and that this change and its amplitude can be detected

with an acceleration sensor.

Methods

Acceleration data synchronized with test stimulation current



Step 1: Data recording setup. A 3 axis accelerometer evaluation board (STEVAL-MKI022V1, ST) housed in a non-conductive printed plastic case (FullCure 830 Vero White, Objet Geometries Ltd) is mounted on the evaluators wrist using a Velcro strap. This sensor is connected to a laptop with in-house developed recording software. This software is also connected to the deep brain stimulation system which provides current during the test stimulations.

Results



Dr. Simone Hemm-Ode T +41 61 4674-796 F +41 61 4674-701 simone.hemm@fhnw.ch

Ashesh Shah T +41 61 4674-413 ashesh.shah@fhnw.ch

Dr. Jerome Coste T +33 4 73 751 002 jcoste@chuclermontferrand.fr





Discussion

- The additional acceleration measurements during the surgery did not increase operation time or the patient's discomfort.
- Sufficient baseline data is necessary for proper identification of acceleration thresholds.
- There is an inherent subjective component in the acceleration analysis because the evaluation is done by the neurologist.
- Further analysis in relation to anatomy could result in better target structures and could raise additional knowledge of the mechanisms of action of DBS

Step 4: The data set with the highest change compared to baseline is identified and the stimulation amplitude corresponding to this data set is defined as the "Quantitatively identified rigidity release" threshold"

Clinical Application

- Clinical study carried out in University Hospital Clermont-Ferrand, France.
- Data was recorded from 9 PD patients who underwent DBS surgery
- A total of 190 test stimulations were performed and data analysed following steps 1 to 4.

Conclusion

- The acceleration of the neurologist's movement is inversely proportional to change in patient's rigidity.
- Acceleration measurements confirm the subjective evaluation, but they seem to be more sensitive (Fig 5).
- Quantitative rigidity evaluation is feasible during DBS surgery.

Acknowledgements

This research has been supported by the Swiss National Science Foundation (SNSF) and the Germaine de Stael program.

Institute for Medical and Analytical Technologies, Head: Prof. Dr.-Ing. Erik Schkommodau, Gründenstrasse 40, 4132 Muttenz, Switzerland. http://www.fhnw.ch/lifesciences/ima